



Oxygen Delivery End Tidal CO₂ Cannula

Ventlab's End Tidal CO₂ Cannula delivers oxygen while capturing an end tidal gas sample, even with simultaneous insufflation of oxygen. The cannula is designed help lessen attenuation of the end tidal sample by oxygen dilution. The End Tidal CO₂ Cannula is non-sterile, latex free.

For detailed product information, visit www.Ventlab.com or contact your Ventlab Sales Representative.

Simultaneously Delivers Oxygen + Captures ETCO₂ Sample



4107
Adult ETCO₂ Cannula / DETAIL

End Tidal CO₂ + Oxygen Delivery Cannula

- Dual nasal cannula simultaneously provides oxygen to both nostrils while obtaining CO₂ sampling during spontaneous breathing
- Designed to prevent mixing of fresh oxygen with CO₂
- 3-Channel tubing prevents kinking
- Includes male/female adapter with luer lock connector
- Available with Fits-All Connector
- Available in three tubing lengths

ETCO₂ + Oxygen Delivery Cannula with Standard Connector



Cat#	Description	Tubing Length	Pkg
4107	Adult ETCO ₂ Cannula	7'	25
4110	Adult ETCO ₂ Cannula	10'	25
4114	Adult ETCO ₂ Cannula	14'	25

ETCO₂ + Oxygen Delivery Cannula with Fits-All Connector



Cat#	Description	Tubing Length	Pkg
4107F	Adult ETCO ₂ Cannula	7'	25
4110F	Adult ETCO ₂ Cannula	10'	25
4114F	Adult ETCO ₂ Cannula	14'	25

