



### SunStim™ Peripheral Nerve Stimulator

- Stimulator output current is fully adjustable from 0-70 mA
- Switch functions:
  - Stand-by position: no stimulus pulses are generated
  - Twitch: 1 pulse per second
  - Tetanus: 100 Hz (100 pulses/s)
  - Built-in train of four: 4 pulses in a period of 2 seconds
- Battery indicator light flashes when battery needs to be replaced
- LED (green) battery light indicates that the unit is on, flashes on low battery
- LED (red) flashes with audible indicator sounds each time a stimulus pulse is generated
- Lightweight: 0.25lbs including battery
- Standard battery: 9 volt alkaline

Cat#	Description	Box
8-1053-60	SunStim™ Peripheral Nerve Stimulator	1

#### Accessories

Cat#	Description	Box
8-1053-76	Ball Electrode and Hollow Pin	1
8-1051-78	Lead 6' Snap and Shielded Pin	1
8-1051-79	Lead 6' Clip and Shielded Pin	1
8-1053-81	SunStim™ Overlay	1

### SunStim™ Plus Peripheral Nerve Stimulator

- Stimulator output current is fully adjustable from 0-70 mA
- Switch functions:
  - Double burst: 2 60-ms bursts of 50 Hz separated by .75 s
  - Twitch: 1 pulse per second
  - Tetanus: 50 Hz (50 pulses/s) or 100 Hz (100 pulses/s)
  - Built-in train of four: 4 pulses in a period of 2 seconds
- Stimulus pulse: square wave monophasic pulses (200- $\mu$ s duration)
- LED (green) battery light indicates that the unit is on, flashes on low battery
- LED (red) flashes with audible indicator sounds each time a stimulus pulse is generated
- Lightweight: 0.25lbs including battery
- Standard battery: 9 volt alkaline

Cat#	Description	Box
8-1053-62	SunStim™ Plus Peripheral Nerve Stimulator	1

#### Accessories

Cat#	Description	Box
8-1053-76	Ball Electrode and Hollow Pin	1
8-1051-78	Lead 6' Snap and Shielded Pin	1
8-1051-79	Lead 6' Clip and Shielded Pin	1
8-1053-82	SunStim™ Plus Overlay	1